

Wound Dressing Kits: *Advanced vs. Basic*

When it comes to wound care, Pensar Medical™ offers two exceptional options to meet your specific needs - Advanced and Basic Wound Dressing Kits. These kits are designed to provide efficient and effective Negative Pressure Wound Therapy (NPWT), ensuring optimal healing outcomes. This one-pager will help you understand the key differences between the two kits to help you make the best choice for your patients.



WoundPro®

Healing. Simplified.

The WoundPro® is an advanced NPWT system renowned for its simplicity and effectiveness. It features the most reliable and powerful pump in the industry to promote accelerated wound healing and enhance patient outcomes.



Advanced Dressing Kits:

Elevating Wound Care Excellence

- ✓ **Ideal for Specialized Care:** The go-to choice for complex wound cases, offering superior control and accuracy.
- ✓ **Includes 10 Kits**, each with:
 - + Hydrophobic foam wound dressings
 - + A dual-lumen StingRay™ Suction Bell
 - + Adhesive drapes
 - + Error-proof, quick-connect tubing connectors
 - + A skin prep pad
- ✓ **Dual-Lumen StingRay™ Suction Bell:** Ensures strong adhesion and precision.
- ✓ **Dual-lumen Tubing:** Reads pressure accurately at the wound site, ideal for high drainage wounds.
- ✓ **Error-Proof Quick-Connect Tubing Connectors:** Prevent leaks and enhance infection control.
- ✓ **Enhanced Canister Compatibility:** Works seamlessly with Advanced Collection Canisters for optimal fluid management.



Basic Dressing Kits:

Reliable and Cost-Effective

- ✓ **Ideal for Routine Care:** Simple and cost-effective, suitable for standard wound care.
- ✓ **Includes 10 Kits**, each with:
 - + Hydrophobic foam wound dressings
 - + A single-lumen StingRay™ Suction Bell
 - + Adhesive drape(s)
 - + An easy-to-use hose clamp
 - + A skin prep pad
- ✓ **Single-Lumen StingRay™ Suction Bell:** Flexible, compact, and comfortable, perfect for difficult-to-place areas.
- ✓ **Easy-To-Position Tubing Connectors:** provides flexibility for occasional minor tubing adjustments.
- ✓ **Easy-To-Use Hose Clamp:** Ensures efficient wound care with airtight changes.
- ✓ **Basic Canister Compatibility:** Compatible with Basic Collection Canisters for straightforward fluid collection.

Both Basic and Advanced Dressing Kits come in small, medium, and large sizes, which is based on the size of the foam included. Small kits come with one drape, whereas medium and large kits have two drapes.

Conclusion

Choose the kit that suits your specific needs. Whether you opt for the precision of the Advanced kits or the economical reliability of the Basic kits, Pensar Medical™ guarantees high-quality products for effective healing and patient comfort.

Contact Us

For more information, please visit our website or contact our team of experts. Experience the Pensar Medical™ difference in wound care today!

Pensar Medical, LLC
Phone: +1 (800) 669-4757
Email: hello@pensarmedical.com