

Our Devices Have a 10-year Track Record of Clinical Success

Advancing the Standard of Wound Healing Care

Pensar Medical is a leader in healthcare innovation, specializing in cutting-edge wound care solutions. Our groundbreaking WoundPro™, WoundPro+™, and MicroDoc® systems enhance patient outcomes and simplify healthcare professionals' daily routines.

Proven to Accelerate Healing of:

- Aesthetic procedures
- Burns
- Complex wounds
- Chronic wounds
- Deep abrasions
- Dermatological procedures
- Diabetic foot ulcers
- Orthopedic surgeries
- Palliative care
- Pressure ulcers
- Skin grafts
- Sports injuries
- Surgical incisions
- Traumatic injuries



WoundPro™

Healing. Simplified.

The WoundPro™ is an advanced NPWT system renowned for its simplicity and effectiveness. It features the most reliable and powerful pump in the industry to promote accelerated wound healing and enhance patient outcomes.

- + **Comprehensive NPWT Therapies:** Offers three NPWT therapy settings – continuous, intermittent, and variable – for versatile patient care.
- + **Industry-Leading Power:** The most powerful pump for consistent results.
- + **Precise Pressure Measurement:** Ensures accurate measurement of the pressure applied at the wound site.



WoundPro+™

Healing. Simplified.

The WoundPro+™ is an advancement from our flagship WoundPro™ system. This NPWT device is now even more compact and lightweight, weighing only half as much as the WoundPro, while retaining its user-friendly interface for effortless operation.

- + **Streamlined Technology:** Our advanced technology ensures accurate tracking, all while maintaining a simplified operation that requires virtually no training, allowing clinicians to focus on patient care.
- + **Proven Reliability:** The WoundPro+™ is powered by the most trusted pump in the market, consistently delivering reliable results.
- + **Enhanced Portability:** With its compact and lightweight design, it's highly portable and ideal for use in hospitals, home health, and long-term care settings.



MicroDoc®

Healing in your pocket

The MicroDoc® is a cutting-edge single-use NPWT system that boasts a non-restrictive, discreet design for easy wearability, all while delivering highly effective wound therapy.

- + **Single-Use & Disposable:** Convenient, hygienic, and maintenance-free.
- + **Effective Wound Healing:** Designed to reduce scarring and infections.
- + **User-Friendly:** Simple one-button operation.
- + **Mobile Patient-Friendly:** Made for active patients, maintaining mobility.
- + **Compact & Discreet:** Tiny, lightweight, wearable on a belt or under clothing.
- + **Whisper-Quiet:** Virtually unnoticed with a gentle, low hum during operation.
- + **Pressure Options:** Three settings (-50, -80, -125 mmHg).



Wound Care Supplies

For Comprehensive Care

Pensar Medical goes beyond standard wound care with our advanced kits, meticulously developed to enhance the efficiency and effectiveness of NPWT.

- + **Collection Canisters:** We offer a variety of collection canister sizes and styles for flexibility in wound care delivery. With a focus on ease of use and disposal, our canisters ensure secure and efficient wound exudate collection.
- + **Foam Dressing Kits:** Our Series Foam Dressing Kits adhere to the highest healthcare standards, offering unparalleled quality and the best materials available to create a secure, sterile environment that promotes effective wound healing while simplifying the dressing change process and minimizing patient care disruptions.

Customer Support

Pensar Medical prioritizes customer support and assistance, ensuring that our valued partners have access to expert guidance when needed.

- ★ Same day order fulfillment ★
- ★ 12-month warranty included ★
- ★ Live support clinical helpline ★

Get In Touch

Ready to elevate your healthcare solutions with Pensar Medical?

www.pensarmedical.com | hello@pensarmedical.com